




PRP TUBES



 407-412, 4th Floor, Span Centre,
South Avenue, Santacruz (W),
Mumbai - 400 054, India.

 +91 22 62396666 / +91 9321759250

 biox@biox.co.in

 www.biox.co.in



CLASSIC PRP TUBE

We would like to introduce ourselves as the marketing partners of Manson PRP tubes.

Manson PRP tubes are used by doctors for Platelet Rich Plasma for Orthopedic, Skin Care, hair regeneration and chronic wound repair.

Being sterile and pyrogen-free the product is ideal for these applications. Most of the other PRP tubes/ ACD tubes available in the market used for these applications are not approved for in vivo application. It's only approved for in vitro i.e. outside the human body.

FEATURES:

- ➔ Usage : Orthopedics , Skin Care , Chronic Wounds , Asthetic Surgery , Hair Regeneration.
- ➔ Materials : PET.
- ➔ Tube Size : 15ml:16 x 125mm.
- ➔ Additives : Sodium Citrate + Separation gel & ACD + Separation gel.
- ➔ Results : Approx 3ml.
- ➔ Centrifugation : 3500 rpm for 8 min.
- ➔ Volume of Blood : 10ml
- ➔ Platelet : 85%
- ➔ Recovery Leukocytes : 20-25%
- ➔ Recovery Erythrocytes : <0.5%



BIOTIN

We would like to introduce ourselves as the marketing partners of Manson Biotin Vials.

Manson Biotin Vials are used with PRP tubes for hair regeneration.

These are sterile , pyrogen free 1ml vials.

