



## BLOOD / FLUID / BLANKET WARMING CABINET



📍 407-412, 4<sup>th</sup> Floor, Span Centre,  
South Avenue, Santacruz (W),  
Mumbai - 400 054, India.

☎ +91 22 62396666 / +91 9321759250

✉ biox@biox.co.in

🌐 www.biox.co.in



## FLUID WARMING CABINET

We would like to introduce ourselves as the manufacturers and supplier of Fluid Warmer Cabinets.

### FEATURES :

- ➔ Comprises forced air warming system for temperature in the cabin.
- ➔ Micro controller based servo proportionate system incorporated to maintain pin point accuracy in the chamber.
- ➔ Special platinum PT-100 Sensor with steel capsule.
- ➔ Digital display provided to indicate set and actual temperature.
- ➔ Accommodates 12, 500 ML fluid bottles.
- ➔ Front door provided with glass for easy viewing.
- ➔ Durable CE marked ceramic heater which has incredible long life.
- ➔ Adjustable temperature range of 32°C to 40°C.
- ➔ Adjustable shelf facility.
- ➔ Illumination lamp provided in the cabinet.



### SPECIFICATION:

MODEL NUMBER	BX-FWC-PETITE	BX-FWC-STANDARD	BX-FWC-GRANDE
Total Capacity	30 Liters	71 Liters	143 Liters
Material Of Construction	MS Powder Coated / Stainless steel (Optional)	MS Powder Coated / Stainless steel (Optional)	MS Powder Coated / Stainless steel (Optional)
Dimensions	20" x 18" x 25" (L x D x H)	20" x 18" x 41" (L x D x H)	20" x 18" x 63" (L x D x H)
Capacity of bottles	7 units of 500ml bottles per rack	7 units of 500ml bottles per rack	7 units of 500ml bottles per rack
Number of racks	2	4	6
Adjustable Temperature Range	32°C to 40°C	32°C to 40°C	32°C to 40°C
Temperature Sensor	Platinum PT-100 sensor with steel capsule	Platinum PT-100 sensor with steel capsule	Platinum PT-100 sensor with steel capsule
Power Supply	220 Volts 50Hz	220 Volts 50Hz	220 Volts 50Hz
Weight	30 Kg	35 Kg	57 Kg