





PATIENT IDENTIFICATION BANDS



 407-412, 4th Floor, Span Centre,
South Avenue, Santacruz (W),
Mumbai - 400 054, India.

 +91 22 62396666 / +91 9321759250

 biox@biox.co.in

 www.biox.co.in



PATIENT IDENTIFICATION BANDS

Patient Identification bands is an important way to reduce patient safety risks. Patient misidentification continues to result in medication errors, transfusion errors, testing errors, wrong person procedures, and the discharge of infants to the wrong families. Today, with as many as one million people harmed by errors during treatment at hospitals each year, identification of patient using wristbands helps minimize medical errors and administrative issues can be handled efficiently.

Bio-X offers optimal and secure patient identification solutions for healthcare facilities.


Colour coding of wristbands facilitates rapid visual recognition of specific issues and can be modified as per the requirement of the organization.

We offer Patient Identification solutions for admissions, blood, mother and baby, paediatrics, disaster response, outpatient and emergency room purposes.

ADULT INSERT CARD BANDS

- Light weight strong, soft vinyl
- Clear shield provides easy, accurate bar code scanning, protection from water & other abrasives
- Non transferable, secure plastic snap provides maximum security
- Information Area: 87mm X 20mm
- Customised printing on request
- Packing: 250 pcs/box


Colour Available: Red, Green, Blue, Yellow, White, Pink, Orange



PEDIATRIC WRITE ON BANDS

- Write on Type Bands
- Soft Vinyl
- Snap Closure
- Information Area: 54mm X 18mm
- Packing: 250 pcs / box


Colour Available: Pink, Blue



MOTHER-INFANT INSERT CARD BANDS

- Mother Baby Insert card type with identification numbers to co-relate the mother with the baby
- Insert Card
- Soft Vinyl
- Snap Closure
- Customised printing on request
- Information Area: 65mm X 21mm (Adult) / 53mm X 15mm (Infant)
- Packing: 250 sets/ box

Colour Available: Pink, Blue, White



PAPER WRITE ON BANDS

- Light weight strong, soft compound material
- Non transferable, secure plastic snap provides maximum security
- Information Area: 87mm X 20mm
- Customised printing on request
- Packing: 250 pcs/box
- Fast and easy option where the patient data can be printed or written with any ink on the bands itself

Colour Available: Light Blue, Dark Blue, Light Pink, Dark Pink, Green, Black, White, Grey, Yellow, Red, Purple, Orange

