




## PATIENT TRANSFER SYSTEM



 407-412, 4<sup>th</sup> Floor, Span Centre,  
South Avenue, Santacruz (W),  
Mumbai - 400 054, India.

 +91 22 62396666 / +91 9321759250

 biox@biox.co.in

 www.biox.co.in




## PATIENT TRANSFER SYSTEM

We would like to introduce ourselves as suppliers of patient transfer system- Easy Move.

This concept of patient transfer is new to India, but it is very popular in the US , UK and Europe for the last 3 decades where lifting patients for transferring is not recommended because it causes the nurses and attendants to suffer back pain/back injuries. The has been the protocol followed in all the emergency departments around the world including Operation Theatres , wards , ICU'S e.t.c : basically where ever patient transfer is done between the OT table to stretcher or stretcher to hospital bed.

Further more it is recommended to use Easy Move because the Patient Transfer take place without lifting the patient's body and especially useful in the case of accident and emergency cases



## FEATURES OF EASY MOVE

Easy Move is made up of a special polycarbonate are plastic which does not break. The length of the board is approximately 170cms and the width is 50cms. The board can be bent over 120°and it will not

### FEATURES :

- ➔ Total Support For The Patient's Body From Head To Toe.
  - ➔ Only 1 Nurse / Attendent Required For A Successful Patient Transfer.
  - ➔ No Risk Of Back Injury & Back Pain To The Nursing Staff Who Aid In Patient Transfer.
  - ➔ Easy Move Lightweight , Radiolucent & Impact Resistant.
  - ➔ Easy Move Made Up Of Special Plastic That Is Unbreakable.
  - ➔ Weight Less Than 4 Kgs , And Can Transfer People With Weight Of Up To 225Kgs.
  - ➔ Life Ruseable If Maintained Properly.
  - ➔ No Change Of Infection To Patient Or Nurse As It Can Be Cleaned With Any Antiseptic Solution Or Mild Disinfectant After Every Procedure If Required.
  - ➔ Help To Bridge The Gap & Height Difference Between Trolley And Hospital Bed. Worldwide Accepted Compared To The Conventional Rollers Which Hurt The Patient And Has A Shelf Life.
- Size : 170 x 50 CM.  
 Weight : 3.35 Kgs.  
 Packing : 1 Pc/Bag.  
 Colour : Blue.