



PATIENT WARMER



PATIENT WARMER

Inadvertent Hypothermia (IHP) is considered to be a grave adverse event for the perioperative patients. While causing the patient discomfort, hypothermia can contribute to complication risks, including morbid cardiac events, surgical site infection, and slow recovering, and eventually results in a prolonged hospital stay.

We would like to introduce ourselves as the channel partners of Keewell Medical technologies revolutionary patient warming system



KPW4000 FORCED-AIR PATIENT WARMING UNIT is a tailor-made & cost-effective solution to prevent peri-operative Inadvertent Hypothermia (IHP) and reduce discomfort caused by cold.

FEATURES :

- ➔ Heating mode : Inflatable convection heating.
- ➔ Display parameters : Set temperature , real-time temperature , heating time , wind speed and fault information.
- ➔ Safety features : Power on self-test, fault category display, over temperature alarm and automatic cut-off of heating and air supply.
- ➔ Disposable blanket : Special double-layer non-woven material, less allergic reaction, antipuncture, anti-crack, water-proof, aseptic packaging.
- ➔ Blower : Built in high-efficiency air filter.
- ➔ Noise : Low.

SPECIFICATION:

Model	:	KPW4000
Input Power	:	a.c.220V/230V, 50Hz/60Hz, 1400VA
Temperature setting value	:	RT (room temperature), 33 °C, 36 °C, 38 °C, 41 °C, 43 °C
Temperature accuracy	:	±1°C
Over temperature power-off protection	:	45°CSoftware/ hardware dual independent protections
Preheating time	:	≤ 2min
Operating mode	:	Continuous
Heating power rated	:	1200W
Type of protection against electric shock	:	Class I
Application part	:	Anti defibrillation BF application part
Moisture protection level	:	IPX2
Dimension	:	(L*W*H) 495*335*465mm
Weight	:	7.5 kg

PATIENT WARMER

Forced-air warming blankets (FAW) are medical devices used to maintain a patient's body temperature during surgery and other medical procedures. These blankets are designed for single use to ensure sterility and prevent cross-contamination between patients.

KPW400 Forced-Air Warming Blankets for Single Use

Full Body KWB 011		Lower Body KWB 012		Upper Body KWB 013	
Upper Arm Extension KWB 014		Lithotomy Underbody KWB 020		Small Full Body KWB 0111	
Small Lower Body KWB 0121		Abdominal Surgery KWB 0151		Full Access Underbody KWB 0171	
Full Access Underbody Upper Arm Extension KWB 0172		Pediatric Underbody KWB 0181		Pediatric Full Access Underbody KWB 0191	





📍 407-412, 4th Floor, Span Centre,
South Avenue, Santacruz (W),
Mumbai - 400 054, India.

☎ +91 22 62396666 / +91 9321759250

✉ biox@biox.co.in

🌐 www.biox.co.in

