




## VEIN VIEWER



 407-412, 4<sup>th</sup> Floor, Span Centre,  
South Avenue, Santacruz (W),  
Mumbai - 400 054, India.

 +91 22 62396666 / +91 9321759250

 [biox@biox.co.in](mailto:biox@biox.co.in)

 [www.biox.co.in](http://www.biox.co.in)



# VEIN VIEWER

We would like to introduce ourselves as the suppliers of Vein Viewer.

## FEATURES:

- ➔ 1.3" DLP display , simple to use.
- ➔ 350g light weight and portable design.
- ➔ 3 different image size are selectable suitable for adult pediatric and neonate.
- ➔ 3 types of display mode , basic , inverse and enhance , meet all the request for a best vasculature map.
- ➔ 10 kinds of color display , design according to different user preference
- ➔ 2 kinds of freehand solution is optional , mobile stand and fixed support.


## SPECIFICATION:

Display	DLP Display
Infrared Wavelength	850nm
Image Resolution	856*480 pixel
Visible Vein Size	≥ 1mm
Penetration Depth	≤8mm
Optimal Imaging Distance	15-25cm
Operating Mode	Basic Mode, Inverse Mode, Enhance Mode
Image Color	Single Color and Mixed Color
Image Size	Large(L), Medium(M), Small(S)
Brightness	High(H), Medium(M), Low(L)
Power Supply	Rechargeable Lithium Battery (3000mAh), DC 5V 2A
Battery Life	2 Hours after fully charged
Net Weight	350g
Dimension	21.2x6.6x6.5cm



### Accuracy Imaging

Advanced software calculations ensure clear image quality and each machine is precisely calibrated at the factory.



### Three Operating Modes

Three operating modes, Basic Mode, Enhance Mode and Inverse Mode, meet the best requirements of all vascular system maps.



**Basic Mode**  
Basic monochromatic mode, clear image



**Enhance Mode**  
Two-color mixing, enhanced appearance



**Inverse Mode**  
Inverse color image, protruding vein

### Three Display sizes

Three different image sizes can be selected for adults, children and neonate.

